

ABL2 protein (STJ505824)
STJ505824**GENERAL INFORMATION**

Product Type	Proteins
Short Description	ABL2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27
Gene Symbol	ABL2
Uniprot ID	ABL2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 1100-1150 on human Tyrosine-protein kinase ABL2/ARG.
Immunogen Region	1100-1150
Specificity	
Immunogen Sequence	