

MLKL protein (STJ505828)

STJ505828

GENERAL INFORMATION

Product Type	Proteins
Short Description	MLKL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74568
Gene Symbol	Mlkl
Uniprot ID	MLKL_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 400-450 on mouse Mixed lineage kinase domain-like protein MLKL.
Immunogen Region	400-450
Specificity	
Immunogen	
Sequence	