

MLKL protein (Phospho) (STJ505830)

STJ505830

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-MLKL-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	197259
Gene Symbol	MLKL
Uniprot ID	MLKL_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 350-400 on human Mixed lineage kinase domain-like protein MLKL
Immunogen Region	350-400
Specificity	
Immunogen	
Sequence	