

EB2 protein (STJ505833)
STJ505833**GENERAL INFORMATION**

Product Type	Proteins
Short Description	EB2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids within region 1-50 on zebra fish Microtubule-associated protein RP/EB family member protein.
Immunogen Region	1-50
Specificity	
Immunogen	
Sequence	