

CGRP protein (STJ505838)
STJ505838**GENERAL INFORMATION**

Product Type	Proteins
Short Description	CGRP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12310
Gene Symbol	Calca
Uniprot ID	CALCA_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 350-410 mouse Calcitonin gene-related peptide 1.
Immunogen Region	350-410
Specificity	
Immunogen Sequence	