

Profilin protein (STJ505842)

STJ505842

GENERAL INFORMATION

Product Type	Proteins
Short Description	Profilin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5216
Gene Symbol	PFN1
Uniprot ID	PROF1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 50-100 on human Profilin 1.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	