

PARD6B protein (STJ505848)

STJ505848

GENERAL INFORMATION

Product Type	Proteins
Short Description	PARD6B-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84612
Gene Symbol	PARD6B
Uniprot ID	PAR6B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	