

## NHERF-1 protein (STJ505849)

STJ505849

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	NHERF-1-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	9368
<b>Gene Symbol</b>	SLC9A3R1
<b>Uniprot ID</b>	NHRF1_HUMAN
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 255-305 on human (+) /H (+) exchange regulatory cofactor NHE-RF1.
<b>Immunogen Region</b>	255-305
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	