

TREML1 protein (STJ505851)

STJ505851

GENERAL INFORMATION

Product Type	Proteins
Short Description	TREML1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	340205
Gene Symbol	TREML1
Uniprot ID	TRML1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 150-200 on human Triggering Receptor Expressed on Myeloid cells-like 1.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	