

ASXL2 protein (STJ505859)

STJ505859

GENERAL INFORMATION

Product Type	Proteins
Short Description	ASXL2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	75302
Gene Symbol	Asxl2
Uniprot ID	ASXL2_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 1200-1275 on mouse Putative Polycomb group protein ASXL2.
Immunogen Region	1200-1275
Specificity	
Immunogen Sequence	