

BMP2 protein (STJ505861)

STJ505861

GENERAL INFORMATION

Product Type	Proteins
Short Description	BMP2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	650
Gene Symbol	BMP2
Uniprot ID	BMP2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 340-396 on human Bone morphogenetic Protein 2.
Immunogen Region	340-396
Specificity	
Immunogen Sequence	