

PARD6B protein (STJ505862)

STJ505862

GENERAL INFORMATION

Product Type	Proteins
Short Description	PARD6B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84612
Gene Symbol	PARD6B
Uniprot ID	PAR6B_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 340-372 on human Partitioning Defective 6 homolog beta.
Immunogen Region	340-372
Specificity	
Immunogen Sequence	