

C1ORF115 protein (STJ505864)

STJ505864

GENERAL INFORMATION

Product Type	Proteins
Short Description	C1ORF115-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution Range	WB: 1:500-1:5, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	226777
Gene Symbol	NA
Uniprot ID	CA115_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Chromosome 1 Open Reading Frame 115.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	