

RHBDF2 protein (STJ505868)

STJ505868

GENERAL INFORMATION

Product Type	Proteins
Short Description	RHBDF2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79651
Gene Symbol	RHBDF2
Uniprot ID	RHBDF2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-50 on human Inactive rhomboid protein 2.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	