

RHBDF2 protein (STJ505871)

STJ505871

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | RHBDF2-protein if for use in research applications. |
| Applications | WB |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel. |
| Conjugation | |
| Purification | |
| Dilution | |
| Range | |
| Formulation | |
| Isotype | |
| Storage | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 79651 |
| Gene Symbol | RHBDF2 |
| Uniprot ID | RHDF2_HUMAN |
| Immunogen | Synthetic peptide taken within amino acid region 1-50 on human Inactive rhomboid protein 2. |
| Immunogen Region | 1-50 |
| Specificity | |
| Immunogen Sequence | |