

GAPDH protein (STJ505877)

STJ505877

GENERAL INFORMATION

Product Type	Proteins
Short Description	GAPDH-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	14433
Gene Symbol	Gapdh
Uniprot ID	G3P_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 150-200 on mouse Glyceraldehyde-3-phosphate dehydrogenase protein.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	