

PDE11A protein (Phospho) (S163) (STJ505879)

STJ505879

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-PDE11A-S163 protein was developed from and has a target region of S163. For use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	140928
Gene Symbol	Pde11a
Uniprot ID	PDE11_RAT
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	