

C1ORF115 protein (STJ505880)

STJ505880

GENERAL INFORMATION

Product Type	Proteins
Short Description	C1ORF115-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 µl per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	79762
Gene Symbol	C1orf115
Uniprot ID	CA115_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 75-125 on human Chromosome 1 Open Reading Frame 115.
Immunogen Region	75-125
Specificity	
Immunogen Sequence	