

TESK2 protein (STJ505889)

STJ505889

GENERAL INFORMATION

Product Type	Proteins
Short Description	TESK2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10420
Gene Symbol	TESK2
Uniprot ID	TESK2_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein.
Immunogen Region	200-250
Specificity	
Immunogen Sequence	