

Androgen Receptor protein (STJ505893)

STJ505893

GENERAL INFORMATION

Product Type	Proteins
Short Description	Androgen Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	367
Gene Symbol	AR
Uniprot ID	ANDR_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 550-600 on Human Androgen Receptor.
Immunogen Region	550-600
Specificity	
Immunogen Sequence	