

CXCR2 protein (STJ505897)

STJ505897

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12765
Gene Symbol	Cxcr2
Uniprot ID	CXCR2_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids within region 445-538 on mouse C-X-C Chemokine Receptot type 2.
Immunogen Region	445-538
Specificity	
Immunogen Sequence	