

CXCR2 protein (STJ505898)

STJ505898

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 Mu l per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12765
Gene Symbol	Cxcr2
Uniprot ID	CXCR2_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	