

CXCR3 protein (STJ505899)

STJ505899

GENERAL INFORMATION

Product Type	Proteins
Short Description	CXCR3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2833
Gene Symbol	CXCR3
Uniprot ID	CXCR3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 175-225 human CXCR3 protein.
Immunogen Region	175-225
Specificity	
Immunogen Sequence	