

NEDD8 protein (STJ505900)

STJ505900

GENERAL INFORMATION

Product Type	Proteins
Short Description	NEDD8-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4738
Gene Symbol	NEDD8
Uniprot ID	NEDD8_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 50-81 on human NEDD8 protein.
Immunogen Region	50-81
Specificity	
Immunogen Sequence	