

SPOP protein (STJ505902)

STJ505902

GENERAL INFORMATION

Product Type	Proteins
Short Description	SPOP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8405
Gene Symbol	SPOP
Uniprot ID	SPOP_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 325-374 on human SPOP protein.
Immunogen Region	325-374
Specificity	
Immunogen Sequence	