

IL18 protein (STJ505906)

STJ505906

GENERAL INFORMATION

Product Type	Proteins
Short Description	IL18-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3606
Gene Symbol	IL18
Uniprot ID	IL18_HUMAN
Immunogen	Recombinant protein corresponding to amino acids within region 37-193 on human IL18 protein.
Immunogen Region	37-193
Specificity	
Immunogen Sequence	