

SLC26A3 protein (STJ505913)

STJ505913

GENERAL INFORMATION

Product Type	Proteins
Short Description	SLC26A3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1811
Gene Symbol	SLC26A3
Uniprot ID	S26A3_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids within region 650-700 on human SCL26A3 protein.
Immunogen Region	650-700
Specificity	
Immunogen Sequence	