

## SLC26A6 protein (STJ505914)

STJ505914

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	SLC26A6-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">65010</a>
<b>Gene Symbol</b>	<a href="#">SLC26A6</a>
<b>Uniprot ID</b>	<a href="#">S26A6_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids within region 200-250 on rat 5-hydroxytryptamine receptor 1A protein.
<b>Immunogen Region</b>	200-250
<b>Specificity</b>	
<b>Immunogen Sequence</b>	