

Mouse Albumin protein (STJ505917)

STJ505917

GENERAL INFORMATION

Product Type	Proteins
Short Description	Mouse Albumin-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	The concentration of the sample is empirically set so that 20 μ l per lane in Western and blotted with corresponding affinity purified antibody will give a clean band with a high signal to noise ratio on the gel.
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11657
Gene Symbol	Alb
Uniprot ID	ALBU_MOUSE
Immunogen	Full length mouse albumin protein.
Immunogen	Full length
Region	
Specificity	
Immunogen	
Sequence	