

Anti-TSLP antibody (50-100) {BIOTIN} (STJ505921)

STJ505921

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TSLP (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85480
Gene Symbol	TSLP
Uniprot ID	TSLP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human Thymic Stromal Lymphopoietin.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	