

Anti-JTB antibody (100-146) {BIOTIN} (STJ505927)

STJ505927

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-JTB (100-146) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10899
Gene Symbol	JTB
Uniprot ID	JTB_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 100-146 on human Jumping Translocation Breakpoint protein.
Immunogen Region	100-146
Specificity	
Immunogen Sequence	