

## Anti-P70 S6 Kinase beta antibody (400-450) {BIOTIN} (STJ505930)

STJ505930

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-P70 S6 Kinase beta (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200 IHC: 1:50-1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">6198</a>
<b>Gene Symbol</b>	<a href="#">RPS6KB1</a>
<b>Uniprot ID</b>	<a href="#">KS6B1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 400-450 on human Ribosomal protein S6 kinase beta protein.
<b>Immunogen Region</b>	400-450
<b>Specificity</b>	
<b>Immunogen Sequence</b>	