

Anti-IFI6 antibody (80-130) (STJ505935)

STJ505935

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IFI6 (80-130) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2537
Gene Symbol	IFI6
Uniprot ID	IFI6_HUMAN
Immunogen	Synthetic peptide take within amino acid region 80-130 on human Interferon alpha-inducible protein 6.
Immunogen	80-130
Region	
Specificity	
Immunogen	
Sequence	