

Anti-SRPK1 antibody (405-455) {BIOTIN} (STJ505942)

STJ505942

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SRPK1 (405-455) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6732
Gene Symbol	SRPK1
Uniprot ID	SRPK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 405-455 on human SRSF protein kinase 1.
Immunogen Region	405-455
Specificity	
Immunogen Sequence	