

Anti-PP2A-alpha antibody (410-486) {BIOTIN} (STJ505945)

STJ505945

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP2A-alpha (410-486) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5525
Gene Symbol	PPP2R5A
Uniprot ID	2A5A_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 410-486 on human Serine/Threonine-Protein Phosphatase 2A, B56-Alpha subunit.
Immunogen Region	410-486
Specificity	
Immunogen Sequence	