

Anti-PP2A-beta antibody (450-497) {FITC} (STJ505949)

STJ505949

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP2A-beta (450-497) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Monkey, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200 IHC: 1:50-1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5526
Gene Symbol	PPP2R5B
Uniprot ID	2A5B_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-497 on human Serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Beta subunit.
Immunogen Region	450-497
Specificity	
Immunogen Sequence	