

Anti-PP2A-epsilon antibody (400-467) {BIOTIN} (STJ505954)

STJ505954

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PP2A-epsilon (400-467) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5529
Gene Symbol	PPP2R5E
Uniprot ID	2A5E_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-467 on human serine/Threonine-Protein Phosphatase 2A 56 kDa regulatory Epsilon subunit.
Immunogen Region	400-467
Specificity	
Immunogen Sequence	