

Anti-CRALBP antibody (260-317) {BIOTIN} (STJ505978)

STJ505978

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CRALBP (260-317) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200 IHC: 1:50-1:250
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6017
Gene Symbol	RLBP1
Uniprot ID	RLBP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 260-317 on human Retinaldehyde-binding protein 1/CRALBP/RLBP1 protein.
Immunogen Region	260-317
Specificity	
Immunogen Sequence	