

## Anti-OTOP1 antibody (370-420) {BIOTIN} (STJ505981)

STJ505981

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-OTOP1 (370-420) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">133060</a>
<b>Gene Symbol</b>	<a href="#">OTOP1</a>
<b>Uniprot ID</b>	<a href="#">OTOP1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 370-420 on human Otopetrin-1.
<b>Immunogen Region</b>	370-420
<b>Specificity</b>	
<b>Immunogen Sequence</b>	