

## Anti-Wilms Tumor Protein antibody {BIOTIN} (STJ505993)

STJ505993

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Wilms Tumor Protein is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications. |
| <b>Applications</b>      | ELISA, IP, WB   |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human, Mouse, Rat   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.5-1.5 µg/µl   |
| <b>Conjugation</b>         | BIOTIN  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:500-1:5, 000<br>ELISA: 1:50, 000<br>IP: 1:200             |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                           |
|---------------------------|---------------------------|
| <b>Gene ID</b>            | <a href="#">7490</a>      |
| <b>Gene Symbol</b>        | <a href="#">WT1</a>       |
| <b>Uniprot ID</b>         | <a href="#">WT1_HUMAN</a> |
| <b>Immunogen</b>          | Wilms tumor protein WT1   |
| <b>Immunogen Region</b>   |                           |
| <b>Specificity</b>        |                           |
| <b>Immunogen Sequence</b> |                           |