

## Anti-PDE3B antibody (260-310) {BIOTIN} (STJ506002)

STJ506002

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-PDE3B (260-310) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Chicken, Human, Monkey, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1.0-1.25 µg/µl
<b>Conjugation</b>	BIOTIN
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:10, 000
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	5140
<b>Gene Symbol</b>	PDE3B
<b>Uniprot ID</b>	PDE3B_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 260-310 on human PDE3B protein. Selective only for PDE3B variant.
<b>Immunogen Region</b>	260-310
<b>Specificity</b>	
<b>Immunogen Sequence</b>	