

Anti-ZAP70 antibody (350-410) {BIOTIN} (STJ506008)

STJ506008

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ZAP70 (350-410) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7535
Gene Symbol	ZAP70
Uniprot ID	ZAP70_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 350-410 on human Tyrosine-protein kinase ZAP-70.
Immunogen Region	350-410
Specificity	
Immunogen Sequence	