

Anti-ASXL3 antibody (2100-2150) {FITC} (STJ506012)

STJ506012

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASXL3 (2100-2150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	211961
Gene Symbol	Asxl3
Uniprot ID	ASXL3_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 2100-2150 on human Putative Polycomb group protein ASXL3.
Immunogen Region	2100-2150
Specificity	
Immunogen Sequence	