

Anti-Phospho-PBXIP1 antibody (100-150) {BIOTIN} (STJ506021)

STJ506021

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-PBXIP1 (100-150) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	229534
Gene Symbol	Pbxip1
Uniprot ID	PBIP1_MOUSE
Immunogen	Phosphorylated synthetic peptide taken within amino acid region 100-150 on mouse Pre-B-cell leukemia transcription factor-interacting protein 1.
Immunogen	100-150
Region	
Specificity	
Immunogen	
Sequence	