

Anti-ADA antibody (50-100) {BIOTIN} (STJ506024)

STJ506024

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADA (50-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100
Gene Symbol	ADA
Uniprot ID	ADA_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 50-100 on human Adenosine Deaminase/Adenosine Aminohydrolase.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	