

Anti-Serine/threonine-protein kinase 4/MST-1 antibody (400-450) (STJ506026)

STJ506026

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine-protein kinase 4/MST-1 (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6789
Gene Symbol	STK4
Uniprot ID	STK4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 400-450 on human Mammalian STE20-like protein kinase 1/STE20-like kinase MST1/STK4.
Immunogen Region	400-450
Specificity	
Immunogen Sequence	