

## Anti-STK3/MST-2 antibody (400-450) (STJ506029)

STJ506029

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-STK3/MST-2 (400-450) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
<b>Applications</b>	ELISA, IP, WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.5-1.5 µg/µl
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB: 1:500 ELISA: 1:50,000 IP: 1:200
<b>Formulation</b>	
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	56274
<b>Gene Symbol</b>	Stk3
<b>Uniprot ID</b>	STK3_MOUSE
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 400-450 on mouse Mammalian STE20-like protein kinase 2/STE20-like kinase MST2/STK3.
<b>Immunogen Region</b>	400-450
<b>Specificity</b>	
<b>Immunogen Sequence</b>	