

Anti-PI3R6 antibody (590-640) (STJ506050)

STJ506050

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PI3R6 (590-640) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	146850
Gene Symbol	PIK3R6
Uniprot ID	PI3R6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 590-640 on human Phosphoinositide 3-kinase regulatory subunit 6.
Immunogen Region	590-640
Specificity	
Immunogen Sequence	