

Anti-EHMT1/GLP1 antibody (950-100) {BIOTIN} (STJ506054)

STJ506054

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EHMT1/GLP1 (950-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	BIOTIN
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50,000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79813
Gene Symbol	EHMT1
Uniprot ID	EHMT1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 950-100 on human Histone-lysine N-methyltransferase EHMT1/GLP1.
Immunogen Region	950-100
Specificity	
Immunogen Sequence	