

Anti-Somatostatin Receptor 2 antibody (300-350) (STJ506056)

STJ506056

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Somatostatin Receptor 2 (300-350) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6752
Gene Symbol	SSTR2
Uniprot ID	SSR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Somatostatin receptor type 2/SSTR2.
Immunogen	300-350
Region	
Specificity	
Immunogen	
Sequence	